

Stops bleeding FAST

CELOXTM

Veterinary



Traumatic injury

Celox, an exceptional haemostat, is finally available to Veterinarians! Extremely effective, safe and easy to use, Celox reliably stops bleeding fast when tackling any trauma bleed situation. You can rest assured, you'll always be prepared when you carry Celox.

Celox is available in two sizes of gauze or granules (just pour, pack and apply pressure). When mixed with blood, Celox forms a gel-like clot in just 30 seconds, independent of the body's normal clotting processes. Made from Chitosan, Celox is safe and naturally metabolised. Unlike some other clotting agents it is non-exothermic so will not burn the patient or care giver. Celox was developed to treat military battlefield wounds and is proven to save lives. Clinical trials by US Marines showed Celox resulted in 100% survival; gave a strong MAP (mean arterial pressure) with a robust clot and no re-bleeding. Independent Veterinary Woundcare Specialist, Georgie Hollis BSc (www.intelligentwoundcare.com) feels that there is a huge comfort knowing you have a sachet of Celox in your bag. Every Veterinarian and horse owner should have it to hand!

Are You Prepared...?



Stake injury – Image courtesy of Greg Quinn, Chine House Veterinary Hospital



Sharp debridement



Tumour removal


CELOXTM is suitable for...

- | Haemostasis following sharp debridement
- | Field trauma
- | Arterial or venous bleeding
- | Severe high pressure bleeding
- | Mass removal
- | Dental bleeds

CELOXTM has been shown to...

- | Significantly reduce blood loss**
- | Reliably stop arterial bleeding**
- | Clot blood containing anticoagulants and Heparin
- | Save lives
- | Save management time for bleeding wounds
- | Be simple and effective

How does Celox work?



Pouring the CELOX into the wound prevents blood loss by forming a gel-like plug as the CELOX links to the surface of red blood cells.

Product ordering information

Order Codes:

- 15g granules FG 08830381
- 5 foot gauze FG 08834151
- 8x8 inch gauze FG 08834141




Medtrade Products Limited
 Electra House
 Crewe Business Park, Crewe,
 Cheshire, UK. CW1 6GL
 Tel: (44) 1270 500019
 Fax: (44) 1270 500045
sales@medtradeproducts.co.uk

* Details of these tests are available at: www.celoxmedical.com
 ** Journal Academic Emergency Medicine Jan 08, Pages 74-81

CELOX™ is easy to use

How to use CELOX™ Gauze



How to use CELOX™ Granules

